



Simec AG

Areal Bleiche West Postfach 413 CH - 4800 Zofingen
Telefon +41 62 752 83 08 Fax +41 62 752 83 09
E-Mail: info@simec.ch www.simec.ch

Test report

Number **C 34812**

Object of analysis: Nandrolone decanoate

Customer: Lightbest International Limited

Samples: 1

Batch: -

Receipt: 05.08.2015

Subject: Assay, Identification

Method(s): HPLC

Sample / Analysis:	Method	Status	Result	Specification	complies
<u>Nandrolone decanoate:</u>					
Identification	HPLC	N	complies	complies	yes
Content	HPLC	N	100.7 %	$\geq 98.5 \%$	yes

Method-status: G = GMP A = accredited V = generally validated P = validated on product N = not validated
E = external lab

Remarks:

Signatures: controlled & released:

Staff:

completed:


20.08.2015

page 1 of 1

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 10:39:51
Vial: 9
Probenidentifikation: STD Nandrolone decanoate
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:30:33
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	25,140	1225852	0,000	0,940
Nandrolone decanoate	36,593	129156842	0,000	99,060
Totals		130382694	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 10-38-37 STD Nandrolone decanoate-Rep1.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 11:33:23
Vial: 9
Probenidentifikation: STD Nandrolone decanoate
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:31:32
Injektvolumen: 20



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
	25,147	1231937	0,000	0,939
Nandrolone decanoate	36,597	129928198	0,000	99,061


Totals		131160135	0,000	100,000
--------	--	-----------	-------	---------

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 11-32-04_STD Nandrolone decanolate-Rep2.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 12:26:46
Vial: 9
Probenidentifikation: STD Nandrolone decanolate
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:31:42
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
	25,133	1222393	0,000	0,936
Nandrolone decanoate	36,590	129313789	0,000	99,064
Totals		130536182	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 12-25-31_STD Nandrolone decanolate-Rep3.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 13:20:09
Vial: 9
Probenidentifikation: STD Nandrolone decanolate
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:31:51
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
	25,143	1218188	0,000	0,935
Nandrolone decanoate	36,587	129000992	0,000	99,065
Totals		130219180	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 13-18-55 STD Nandrolone decanolate-Rep4.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 14:13:32
Vial: 9
Probenidentifikation: STD Nandrolone decanoate
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:32:04
Injektvolumen: 20



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
	25,157	1242845	0,000	0,938
Nandrolone decanoate	36,607	131235658	0,000	99,062


Totals		132478503	0,000	100,000
--------	--	-----------	-------	---------

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 14-12-20_STD Nandrolone decanoate-Rep5.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 15:06:56
Vial: 9
Probenidentifikation: STD Nandrolone decanolate
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:36:22
Injektvolumen: 20



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
	25,167	1242161	0,000	0,942
Nandrolone decanoate	36,630	130640455	0,000	99,058


Totals		131882616	0,000	100,000
--------	--	-----------	-------	---------

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 15-05-42_STD Nandrolone decanolate-Rep6.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 16:00:31
Vial: 1
Probenidentifikation: Blank
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:36:32
Injektvolumen: 20



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
Nandrolone decanoate			0,000 BDL	


Totals				
--------	--	--	--	--

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 15-59-11_Blank.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 16:54:00
Vial: 10
Probenidentifikation: C34812 A
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:36:43
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	36,630	151370427	0,000	100,000
Totals		151370427	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 16-52-44_C34812 A-Rep1.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 17:47:27
Vial: 10
Probenidentifikation: C34812 A
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:36:52
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	36,600	152231615	0,000	100,000
Totals		152231615	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 17-46-12_C34812 A-Rep2.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 18:40:56
Vial: 10
Probenidentifikation: C34812 A
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:37:03
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	36,573	151769419	0,000	100,000
Totals		151769419	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 18-39-40_C34812 A-Rep3.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 19:34:22
Vial: 11
Probenidentifikation: C34812 B
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:37:14
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	36,643	123716418	0,000	100,000
Totals		123716418	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 19-33-05_C34812 B-Rep1.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 20:27:49
Vial: 11
Probenidentifikation: C34812 B
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:37:26
Injektvolumen: 20



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	36,590	124222209	0,000	100,000
Totals		124222209	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 20-26-30_C34812 B-Rep2.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
Analysendatum: 19.08.2015 21:21:14
Vial: 11
Probenidentifikation: C34812 B
Operator: admin
Chrom.-Beding. ACN / Wasser, gradient
Temperatur: 40 °C
Säule Nr. 88, 250 x 4,6, C18
Detektion: UV 240nm

Druckdatum: 20.08.2015 08:37:38
Injektvolumen: 20



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	36,650	123643516	0,000	100,000
Totals		123643516	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2015\August 2015\Diverse
Steroide\18,08,15\19.08.2015 21-20-00_C34812 B-Rep3.dat
Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met